

Pittsburgh General - Biomed

Asset Risk Defaults with PM Frequency

Category Code	Category Description	Sub Category Code	Sub Category Description	Risk Strategy	Default Risk #	Frequency
#ECGM	#ECGM					
		ECG M0	non critical monitor	Clinical Equipment Risk Classification	10.67	Annually
		ECG MO	ECG MONITOR	Clinical Equipment Risk Classification	13.00	Annually
AIR PA	Self Contained (SCBA)					
		AIR PA	Emergency air supply for MRI exam rooms	Clinical Equipment Risk Classification	8.67	Monthly
ALARM0	Equipment that is alarmed					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
AMISC	AMISC					
		MISC0	Not a equipment category	Clinical Equipment Risk Classification	3.00	Not Required
AMP	AMP					
		AMPLI0	AMPLIFIER RACK	Clinical Equipment Risk Classification	4.00	Not Required
		AMPLIF	Amplifier for overhead paging system	Clinical Equipment Risk Classification	3.00	Not Required
		PACEM0	Non- Critical Pacemaker device	Clinical Equipment Risk Classification	9.33	Annually
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
ANALYZX	Analyzer					
		AGM NC	GAS ANALYZER	Clinical Equipment Risk Classification	11.67	Annually
		ANALY0	Large Analyzers	Clinical Equipment Risk Classification	10.67	Annually
		ERGON0	exercise device	Clinical Equipment Risk Classification	13.00	Annually
		EVENT	Event Marker	Clinical Equipment Risk Classification	8.33	Annually
		PACEM0	Non- Critical Pacemaker device	Clinical Equipment Risk Classification	9.33	Annually
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
AUDIO	Audiovisual Equipment					
		A/V EQ	A/V EQUIP	Clinical Equipment Risk Classification	6.67	Annually
		AMPLIF	Amplifier for overhead paging system	Clinical Equipment Risk Classification	3.00	Not Required
		CAMERA	Photo device or attachments	Clinical Equipment Risk Classification	7.00	Annually
		MIXER/	audio equipment	Clinical Equipment Risk Classification	8.33	Not Required
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		STEREO	Stereo receiver	Clinical Equipment Risk Classification	4.33	Not Required

Report Requirements:

1. Asset Risk Management must be enabled
2. Asset Sub Categories must have defaults.
3. Risk Sub Category referring to PM needs must have "preventive" or "PM" in the name.
4. Only Asset Sub Category Codes with defaults are listed.

Category Code	Category Description	Sub Category Code	Sub Category Description	Risk Strategy	Default Risk #	Frequency
		TELECA	Endoscope telecam	Clinical Equipment Risk Classification	12.33	Semi-Annually
		VIDEO5	VIDEO PROCESSOR	Clinical Equipment Risk Classification	13.00	Semi-Annually
<u>AUDITO</u>	<u>AUDITORY</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		TEST E	for equip pm troubleshooting/calibration	Clinical Equipment Risk Classification	2.67	Not Required
<u>AV</u>	<u>Audio Visual</u>					
		CONFER	CONFERENCE ROOM	Clinical Equipment Risk Classification	3.00	Not Required
		MISC0	Not a equipment category	Clinical Equipment Risk Classification	3.00	Not Required
		VTR	Video tape recorder - One inch	Clinical Equipment Risk Classification	4.00	Not Required
<u>BOX</u>	<u>Box</u>					
		REINF	REINF BOX	Clinical Equipment Risk Classification	5.67	Not Required
		STENT	STENT MAKER	Clinical Equipment Risk Classification	6.00	Not Required
<u>BRONK</u>	<u>Bronk</u>					
		LECTUR	LECTURE SCOPE	Clinical Equipment Risk Classification	8.00	Not Required
<u>BYPASS</u>	<u>Bypass</u>					
		AGM	AGM	Clinical Equipment Risk Classification	15.00	Annually
		ERC	Electrical Remote Clamp for ECMO	Clinical Equipment Risk Classification	16.33	Semi-Annually
		INCUBA	Specimen processing in lab areas	Clinical Equipment Risk Classification	9.67	Annually
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		WATER2	Non Critical water module	Clinical Equipment Risk Classification	9.67	Not Required
<u>CAMER1</u>	<u>Camera</u>					
		CAMERA	Photo device or attachments	Clinical Equipment Risk Classification	7.00	Annually
<u>CD BUR</u>	<u>CD Burner</u>					
		MISC0	Not a equipment category	Clinical Equipment Risk Classification	3.00	Not Required
<u>CHARGE</u>	<u>Charger</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		VAD	Left Ventricular Assist Device	Clinical Equipment Risk Classification	18.00	Semi-Annually
		VAD NC	Ventricular Assist Device Non Critical	Clinical Equipment Risk Classification	10.33	Annually
<u>COMPROX</u>	<u>Compression Device</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required

Report Requirements:

1. Asset Risk Management must be enabled
2. Asset Sub Categories must have defaults.
3. Risk Sub Category referring to PM needs must have "preventive" or "PM" in the name.
4. Only Asset Sub Category Codes with defaults are listed.

Category Code	Category Description	Sub Category Code	Sub Category Description	Risk Strategy	Default Risk #	Frequency
<u>COMPUTX</u>	<u>Computer</u>					
		DENT	Dental X-ray Equipment	Clinical Equipment Risk Classification	12.67	Annually
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		SCC	SCC	Clinical Equipment Risk Classification	12.67	Annually
<u>COVER0</u>	<u>Cover</u>					
		COVER0	Cover	Clinical Equipment Risk Classification	8.00	Annually
<u>CPAP</u>	<u>CPAP</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>CPU</u>	<u>CPU</u>					
		ARRYTH	ARRYTHMIA CPU	Clinical Equipment Risk Classification	14.33	Semi-Annually
<u>CYCLOT</u>	<u>CN Cyclotron Lab areas</u>					
		ALL-M4	Lab equipment	Clinical Equipment Risk Classification	10.33	Not Required
		MISC0	Not a equipment category	Clinical Equipment Risk Classification	3.00	Not Required
<u>DEFIB</u>	<u>DEFIB</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>DENTALX</u>	<u>DENTAL</u>					
		DENTAL	Dental Lab Engine	Clinical Equipment Risk Classification	12.33	Annually
<u>DIALYS</u>	<u>DIALYSIS</u>					
		CONCEN	Acid and Bicarb Delivery Systems for Dialysis	Clinical Equipment Risk Classification	3.67	Annually
		MISC	Misc. Equipment or System which is Non Critic	Clinical Equipment Risk Classification	3.00	Not Required
		UNIT S	Dialysis Unit Support	Clinical Equipment Risk Classification	3.00	Not Required
<u>ECG</u>	<u>ECG</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>ECGHP</u>	<u>EGC HP</u>					
		AGM	AGM	Clinical Equipment Risk Classification	15.00	Annually
		CPU	CPU	Clinical Equipment Risk Classification	11.33	Annually
<u>ECGPHI</u>	<u>ECG Phil</u>					
		PDS	PDS	Clinical Equipment Risk Classification	11.33	Annually
<u>EEG</u>	<u>EEG</u>					
		EEG EQ	EEG EQUIPMENT	Clinical Equipment Risk Classification	12.00	Annually

Report Requirements:

1. Asset Risk Management must be enabled
2. Asset Sub Categories must have defaults.
3. Risk Sub Category referring to PM needs must have "preventive" or "PM" in the name.
4. Only Asset Sub Category Codes with defaults are listed.

Category Code	Category Description	Sub Category Code	Sub Category Description	Risk Strategy	Default Risk #	Frequency
		EEG RE	EEG READER	Clinical Equipment Risk Classification	11.00	Semi-Annually
		EEG SY	EEG SYSTEM	Clinical Equipment Risk Classification	12.00	Semi-Annually
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		SYSTE0	system power supply for Sleep Lab systems	Clinical Equipment Risk Classification	8.33	Annually
<u>EIM</u>	<u>EIM</u>					
		NONE	not equipment	Clinical Equipment Risk Classification	3.00	Not Required
<u>EIO</u>	<u>EIO</u>					
		CAMERA	Photo device or attachments	Clinical Equipment Risk Classification	7.00	Annually
		MISC0	Not a equipment category	Clinical Equipment Risk Classification	3.00	Not Required
<u>ELECT3X</u>	<u>ELECTROSURGICAL</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>EXERCIX</u>	<u>EXERCISE EQUIPMENT</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		THERA1	THERAPY DEVICE	Clinical Equipment Risk Classification	13.67	Annually
<u>EXR</u>	<u>EXR</u>					
		MISC0	Not a equipment category	Clinical Equipment Risk Classification	3.00	Not Required
<u>HVAC</u>	<u>HVAC</u>					
		CHILLE	CHILLED WATER PUMP	Clinical Equipment Risk Classification	3.33	Not Required
<u>INS</u>	<u>INS</u>					
		LASER0	Surgical or therapeutic class III and IV lasers	Clinical Equipment Risk Classification	16.33	Annually
		LIGHT1	Imaging device that uses particular light	Clinical Equipment Risk Classification	12.33	Annually
		MISC0	Not a equipment category	Clinical Equipment Risk Classification	3.00	Not Required
		MISCEL	worst case for this grouping is a score of 11	Clinical Equipment Risk Classification	11.00	Annually
		NONE	not equipment	Clinical Equipment Risk Classification	3.00	Not Required
		PACS 2	electronic storage and display of images	Clinical Equipment Risk Classification	4.00	Not Required
		STAINE	Lab equipment	Clinical Equipment Risk Classification	8.00	Not Required
<u>IRRIGAX</u>	<u>IRRIGATION</u>					
		IRRIG0	Irrigation Module attached to Electrosurgical U	Clinical Equipment Risk Classification	13.33	Semi-Annually
<u>LAB</u>	<u>LAB</u>					
		PV LAB	Portable Vascular Lab	Clinical Equipment Risk Classification	12.33	Semi-Annually

Report Requirements:

1. Asset Risk Management must be enabled
2. Asset Sub Categories must have defaults.
3. Risk Sub Category referring to PM needs must have "preventive" or "PM" in the name.
4. Only Asset Sub Category Codes with defaults are listed.

Category Code	Category Description	Sub Category Code	Sub Category Description	Risk Strategy	Default Risk #	Frequency
<u>LAB EQ</u>	<u>LAB EQUIPMENT</u>	SYNTHE	Processes radioisotopes coming from the cycl	Clinical Equipment Risk Classification	5.67	Annually
<u>LAMP</u>	<u>LAMP</u>	HEMATO	Used in labs	Clinical Equipment Risk Classification	10.00	Annually
<u>LIGHT1</u>	<u>LIGHT</u>	PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>MACHIN</u>	<u>MACHINE</u>	PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		CONSOL	CONSOLE	Clinical Equipment Risk Classification	14.33	Semi-Annually
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>MECV</u>	<u>MECV</u>	THERA1	THERAPY DEVICE	Clinical Equipment Risk Classification	13.67	Annually
<u>METER</u>	<u>METER</u>	MISC	Misc. Equipment or System which is Non Critic	Clinical Equipment Risk Classification	3.00	Not Required
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		SOUND	SOUND METER	Clinical Equipment Risk Classification	9.00	Annually
<u>MISC</u>	<u>Miscellaneous</u>	SUDORO	SUDOROMETER	Clinical Equipment Risk Classification	11.00	Annually
		MISC	Misc. Equipment or System which is Non Critic	Clinical Equipment Risk Classification	3.00	Not Required
		MISC0	Not a equipment category	Clinical Equipment Risk Classification	3.00	Not Required
<u>MIXER</u>	<u>MIXER</u>	RESPM	RESPM	Clinical Equipment Risk Classification	11.67	Annually
<u>MODEM</u>	<u>MODEM</u>	CONCEN	Acid and Bicarb Delivery Systems for Dialysis	Clinical Equipment Risk Classification	3.67	Annually
<u>MODULE</u>	<u>MODULE</u>	PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>MONITOX</u>	<u>MONITOR</u>	PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		DIRECT	Monitor	Clinical Equipment Risk Classification	13.00	Semi-Annually
		EMG MA	EMG MACHINE	Clinical Equipment Risk Classification	14.33	Annually

Report Requirements:

1. Asset Risk Management must be enabled
2. Asset Sub Categories must have defaults.
3. Risk Sub Category referring to PM needs must have "preventive" or "PM" in the name.
4. Only Asset Sub Category Codes with defaults are listed.

Category Code	Category Description	Sub Category Code	Sub Category Description	Risk Strategy	Default Risk #	Frequency
		MISC0	Not a equipment category	Clinical Equipment Risk Classification	3.00	Not Required
		MONITE	MONITOR-EXERCISE	Clinical Equipment Risk Classification	12.33	Annually
		MONITH	Wavemap Pressure Guide - Pressure/Flow Sy	Clinical Equipment Risk Classification	13.00	Semi-Annually
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		PRESS3	Used to monitor the patient cavity pressure du	Clinical Equipment Risk Classification	9.00	Annually
		TEMPE0	Monitor	Clinical Equipment Risk Classification	6.33	Semi-Annually
<u>MOTOR</u>	<u>MOTOR</u>					
		DERMAT	skin graphing	Clinical Equipment Risk Classification	12.00	Annually
		EYE MA	EYE MAGNET	Clinical Equipment Risk Classification	14.00	Annually
		HAND E	HAND ENGINE	Clinical Equipment Risk Classification	12.00	Annually
<u>NEBULIX</u>	<u>NEBULIZER</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>NETWORK</u>	<u>NETWORK</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>NIBP</u>	<u>NIBP</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>NON EQ</u>	<u>Non Equipment Class</u>					
		ALL-IN	Equip Id code 351	Clinical Equipment Risk Classification	9.67	Annually
		ALL-LE	Equip ID code 352	Clinical Equipment Risk Classification	10.33	Annually
		ALL-OP	Equip ID 350	Clinical Equipment Risk Classification	9.33	Not Required
		DRILL	Machining Tool	Clinical Equipment Risk Classification	4.00	Not Required
<u>NONE</u>	<u>NONE</u>					
		NONE	not equipment	Clinical Equipment Risk Classification	3.00	Not Required
<u>OPT</u>	<u>OPT</u>					
		MISC0	Not a equipment category	Clinical Equipment Risk Classification	3.00	Not Required
		MISCEL	worst case for this grouping is a score of 11	Clinical Equipment Risk Classification	11.00	Annually
<u>PC INA</u>	<u>PC INAC</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>PERIPH</u>	<u>PERIPHERAL</u>					
		INTERF	interface non-critical	Clinical Equipment Risk Classification	11.33	Not Required

Report Requirements:

1. Asset Risk Management must be enabled
2. Asset Sub Categories must have defaults.
3. Risk Sub Category referring to PM needs must have "preventive" or "PM" in the name.
4. Only Asset Sub Category Codes with defaults are listed.

Category Code	Category Description	Sub Category Code	Sub Category Description	Risk Strategy	Default Risk #	Frequency
		KEYBOA	KEYBOARD	Clinical Equipment Risk Classification	3.33	Not Required
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		PDA	Personal Data Assistant	Clinical Equipment Risk Classification	8.67	Not Required
		PRINTE	Paper printer	Clinical Equipment Risk Classification	7.00	Annually
		TERMIN	TERMINAL	Clinical Equipment Risk Classification	4.67	Not Required
<u>POD</u>	<u>POD</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		POD	POD	Clinical Equipment Risk Classification	11.00	Annually
<u>POWER0X</u>	<u>POWER SUPPLY</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		POWERA	Power supply test instrument	Clinical Equipment Risk Classification	3.00	Not Required
<u>PRINTEX</u>	<u>Printer</u>					
		COPIER	Buisness Copier	Clinical Equipment Risk Classification	6.33	Not Required
		PRINTE	Paper printer	Clinical Equipment Risk Classification	7.00	Annually
<u>PROGRAX</u>	<u>PROGRAMMER</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>PULOX</u>	<u>PULOX</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>PUMP</u>	<u>Pump</u>					
		ENDOFL	Used to inflate the patient cavity during surger	Clinical Equipment Risk Classification	9.00	Annually
		HYDRO1	Water jet for tissue removal	Clinical Equipment Risk Classification	9.33	Annually
		LAPARO	Laparoflator	Clinical Equipment Risk Classification	14.67	Annually
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		REPEAT	Repeater pump for pharmacy production IV's	Clinical Equipment Risk Classification	9.33	Annually
		WASTE1	Patient fluid and smoke treatment system	Clinical Equipment Risk Classification	6.67	Annually
<u>RACK</u>	<u>RACK</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>RAD</u>	<u>RAD</u>					
		INS TR	Instrument Tracking System	Clinical Equipment Risk Classification	9.00	Annually
		MISC	Misc. Equipment or System which is Non Critic	Clinical Equipment Risk Classification	3.00	Not Required

Report Requirements:

1. Asset Risk Management must be enabled
2. Asset Sub Categories must have defaults.
3. Risk Sub Category referring to PM needs must have "preventive" or "PM" in the name.
4. Only Asset Sub Category Codes with defaults are listed.

Category Code	Category Description	Sub Category Code	Sub Category Description	Risk Strategy	Default Risk #	Frequency
		MISC0	Not a equipment category	Clinical Equipment Risk Classification	3.00	Not Required
		NONE	not equipment	Clinical Equipment Risk Classification	3.00	Not Required
		RESPIR	Respiratory Gating of Imaging Equip.	Clinical Equipment Risk Classification	9.00	Not Required
		ULTRA1	Table used in a diagnostic ultrasound room	Clinical Equipment Risk Classification	4.67	Not Required
		ULTRAB	Automated Breast Volume Scanner Transducer	Clinical Equipment Risk Classification	11.67	Annually
<u>RAT</u>	<u>RAT</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>RECORDX</u>	<u>RECORDER</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>REFRIO</u>	<u>REFRIGERATION</u>					
		SLUSH	SLUSH BASIN	Clinical Equipment Risk Classification	13.00	Annually
<u>SCALE</u>	<u>SCALE</u>					
		WEIGHT	WEIGHT MACHINE	Clinical Equipment Risk Classification	12.00	Annually
<u>SENSO0</u>	<u>SENSORY</u>					
		SENSOR	SENSORY	Clinical Equipment Risk Classification	13.00	Annually
<u>SENSORX</u>	<u>SENSOR</u>					
		CO2 MO	CO2 MONITOR	Clinical Equipment Risk Classification	15.67	Annually
		HEATER	lab equip	Clinical Equipment Risk Classification	7.00	Not Required
<u>SERVER</u>	<u>SERVER</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>SIMULAX</u>	<u>SIMULATOR</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>STIMU0</u>	<u>STIMULATOR</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
<u>STROBE</u>	<u>STROBE</u>					
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		STROBE	STROBE	Clinical Equipment Risk Classification	11.00	Annually
<u>SURGER</u>	<u>SURGERY MISC</u>					
		INJEC1	Dutch Ophthalmic Research Center System	Clinical Equipment Risk Classification	16.33	Semi-Annually
<u>SWITCH</u>	<u>SWITCH</u>					

Report Requirements:

1. Asset Risk Management must be enabled
2. Asset Sub Categories must have defaults.
3. Risk Sub Category referring to PM needs must have "preventive" or "PM" in the name.
4. Only Asset Sub Category Codes with defaults are listed.

Category Code	Category Description	Sub Category Code	Sub Category Description	Risk Strategy	Default Risk #	Frequency
<u>SYSTEM</u>	<u>SYSTEM</u>	PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		SWITCH	switch box	Clinical Equipment Risk Classification	9.33	Annually
<u>TABLE</u>	<u>TABLE</u>	ACTIVA	Stryker Activation System	Clinical Equipment Risk Classification	6.67	Not Required
		PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		QMT SY	QMT SYSTEM	Clinical Equipment Risk Classification	13.67	Semi-Annually
<u>TELEMEX</u>	<u>TELEMETRY</u>	EXAM T	Patient exam table	Clinical Equipment Risk Classification	7.33	Semi-Annually
		THERA0	Table used in physical therapy unit	Clinical Equipment Risk Classification	11.00	Annually
<u>TEST</u>	<u>TEST</u>	PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		TEST	TEST	Clinical Equipment Risk Classification	1.00	
<u>THERAP</u>	<u>Demistifierr</u>	THERAP	Demistifier	Clinical Equipment Risk Classification	6.00	Annually
<u>TRAINE</u>	<u>TRAINER</u>	THERA1	THERAPY DEVICE	Clinical Equipment Risk Classification	13.67	Annually
<u>TRAINIX</u>	<u>TRAINING</u>	COVER0	N/A	Clinical Equipment Risk Classification	8.00	Annually
<u>ULTRA0</u>	<u>Ultrasound</u>	DENTA0	Cutting bone tissue and gums	Clinical Equipment Risk Classification	16.00	Annually
<u>UNKNOWN</u>	<u>Unknown</u>	PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		COMPUT	COMPUTER CRT	Clinical Equipment Risk Classification	11.33	Annually
		CPU	CPU	Clinical Equipment Risk Classification	11.33	Annually
<u>UPS</u>	<u>UPS</u>	HAND 0	Hand drill	Clinical Equipment Risk Classification	8.00	Annually
<u>VACUUMX</u>	<u>Vacuum</u>	PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		SMOKE2	Evacuator used with electrosurgical generator	Clinical Equipment Risk Classification	6.33	Semi-Annually

Report Requirements:

1. Asset Risk Management must be enabled
2. Asset Sub Categories must have defaults.
3. Risk Sub Category referring to PM needs must have "preventive" or "PM" in the name.
4. Only Asset Sub Category Codes with defaults are listed.

Category Code	Category Description	Sub Category Code	Sub Category Description	Risk Strategy	Default Risk #	Frequency
<u>VALVE</u>	<u>Valve</u>	VACUU7	Vacuum System	Clinical Equipment Risk Classification	7.67	Not Required
		SWITC2	pneumatic switching valve	Clinical Equipment Risk Classification	10.33	Annually
		VALVE0	VALVE	Clinical Equipment Risk Classification	5.67	Not Required
<u>VESTIB</u>	<u>Vestibular</u>	PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		VESTIB	VESTIBULAR EQU	Clinical Equipment Risk Classification	12.67	Annually
<u>VIDEOX</u>	<u>Video</u>	VIDEO1	distribution of video in diagnostic Radiology ar	Clinical Equipment Risk Classification	4.33	Not Required
<u>WARMERX</u>	<u>Warmer</u>	PC-4	Equip removed from Inventory	Clinical Equipment Risk Classification	4.00	Not Required
		POWER6	Battery pack	Clinical Equipment Risk Classification	5.00	Annually
<u>WASHERX</u>	<u>Washer</u>	PLATE2	Lab equipment	Clinical Equipment Risk Classification	8.67	Annually

Report Requirements:

1. Asset Risk Management must be enabled
2. Asset Sub Categories must have defaults.
3. Risk Sub Category referring to PM needs must have "preventive" or "PM" in the name.
4. Only Asset Sub Category Codes with defaults are listed.